## SBG RD168 WTx metadata

In addition to the metadata below, fastqs will be tagged with the library batch and sequencing run number, in this format: "Seq0x\_Lib0x" (eg. Seq01\_Lib01)

	FIXED VALUE	VARIABLE
Experimental strategy	Total RNA-Seq	
Library ID		CGR sample ID_with version # if not v1 (eg. "SC080048_v1")
Platform	Illumina HiSeq	
Platform Unit ID		Runxx_Lanexx (eg "Run01_Lane01")
File Segment Number		Lane number (integer only, eg "1")
Quality Scale	illumina 18	
Paired-end		Set to 1 or 2 (R1 = 1, R2 = 2)
Reference genome	ucsc.hg19	

	FIXED VALUE	VARIABLE
Sample ID		CGR sample ID
Sample type		Tumor Tissue or Normal Tissue
Investigation	-	-
Species	Homo sapiens	
Batch number		Libxx_Seqxx (Library batch # and Sequencing Run #
Case ID		Patient ID (from Nationwide spreadsheet; eg UA2124)

leaving blank for now - may use it later on for variables like radiation dose, etc

Green fields are consistent between all samples